




Centrale di evacuazione gas anestetici GSE

Anaesthetic gas evacuation system model GSE



		Modello/Model			
		GSE 30	GSE 70	GSE 100	GSE 130
Dimensioni (mm) Dimensions (mm)	A	1" G	1¼" G	1½" G	1½" G
	B	1" G	1¼" G	1½" G	1½" G
	C	700	800	900	960
	D	500	550	650	700
	E	35	35	50	50
	F	130	140	150	110
	G	240	260	300	400
	H	600	600	600	700
	L	650	650	650	750
Peso/Weight	kg	60	80	100	110